

ILLUMINOSS® — Transforming Osteoporotic and Pathologic Fracture Repair

*Advancing fracture fixation through
minimally invasive, conforming
intramedullary implants*



ILLUMINOSS®

A NEW TREATMENT PARADIGM

for fracture repair

The IlluminOss® Solution delivers a completely new approach to fracture repair which may promote patients' return to mobility and activities of daily living.

IlluminOss addresses key challenges when treating fractures in poor quality bone by:

- *Enabling a minimally invasive surgical approach.*^{1,2} IlluminOss requires a small incision for access and implant insertion — minimizing potential soft tissue damage.
- *Conforming to the intramedullary canal,*^{1,3,4} the implant expands to fit the unique geometry of the patient's canal.
- *Delivering rapid longitudinal and rotational stability* at the fracture site and across the length of the implant.^{1,2,3,4}
- *Accommodating the use of ancillary fixation devices along the length of the implant.* IlluminOss augments use of conventional plates, screws and nails for enhanced stability in compromised bone.^{1,2,3,4}




¹ Heck S, Gick S, Rabiner R, Penning D, New Strategy in Geriatric Traumatology — First use of an intramedullary photodynamic polymer in the humerus. VSOU Baden-Baden (South German Orthopedic Congress). April 28-May 1, 2012.

² Vegt PA, Muir JM, Block JE. The Photodynamic Bone Stabilization System: a minimally invasive, percutaneous intramedullary polymeric osteosynthesis for simple and complex long bone fractures. Medical Devices: Evidence and Research. 2014; 7:453-461

³ Heck S, Gick S, Penning D, Intramedullary polymer implant with angular stability for minimally invasive repair of pathologic fractures. Meeting of the Associations of Surgeons — Cologne, Germany.



Balloon is infused with liquid monomer causing it to expand and conform to the canal



The monomer is cured on demand by visible light



Implant provides longitudinal and rotational stability along the length of the implant

⁴ Heck S, Gick S, Penning D, Minimally invasive stabilization of upper limb pathological fractures with an intramedullary polymer. AAOS (American Academy of Orthopedic Surgeons). New Orleans, March 11-15, 2014

Indications

The IlluminOss System is CE Mark approved for light to low, load-bearing bones. It is indicated for use in fracture alignment reduction. It provides stabilization for bone fractures using a percutaneous technique in which the bone is not subjected to significant weight bearing forces. This product has the CE mark and is available for sale in the EU. This product is investigational and not for sale in the United States.

Mobility and independence regained

"Treatment with IlluminOss contributed to a faster return to baseline mobility for this 88-year-old patient with a fibula fracture. She was able to resume walking one day post surgically."

~ Thomas Gausepohl, MD, PhD,
Klinikum Wetzlar, Germany



Scan the QR code to see
a fibula fracture procedure
and post-procedure result.



ILLUMINOSS

IlluminOss Medical, Inc.
993 Waterman Avenue
East Providence, Rhode Island
02914 USA
Tel: +1 (401) 714 0008
Fax: +1 (401) 714 0009
info@illuminoss.com
illuminoss.com

IlluminOss Medical BV
Avenue Ceramique 223
6221 KX Maastricht, Netherlands
info@illuminoss.com

©2016 Illuminoss Medical, Inc.

IlluminOss Medical GmbH
Mühlenhof 7-9
40721 Hilden, Germany
info@illuminoss.com

All rights reserved 900478.B