**DESCRIPTION**

IlluminOss offers a variety of reusable surgical instruments to assist the surgeon with the implantation procedure when, in their opinion, other hand tools may not be convenient. These surgical instruments must be sterilized before use. The information provided is for the surgical instruments only, and do not apply to the Bone Stabilization System catheters.

**INDICATIONS FOR USE**

The IlluminOss Photodynamic Bone Stabilization System is indicated for use in fracture alignment reduction. It provides stabilization for bone fractures using a minimally invasive technique in which the bone is not subjected to weight bearing forces. All bones in which the IlluminOss device is intended to treat include: metacarpal, phalanges, clavicle, radius, ulna, distal radius, olecranon and fibula.

**CONTRAINDICATIONS**

- Patients who are considered skeletally immature.
- Presence of active or incompletely treated infections that could involve the site where the device will be implanted.
- Patients allergic to any of the implant materials, or to dental glue.
- Patients whose intramedullary canal at site of fracture measures smaller than 3.5mm in diameter.
- Uncooperative patient or patient with neurologic disorder, incapable of following directions.
- Metabolic disorders which may impair bone formation.
- Osteomalacia.
- Distant foci of infections which may spread to the implant site.
- Marked bone loss or bone resorption.
- Vascular insufficiency, muscular atrophy, or neuro-muscular disease.
- Open fractures with severe contamination.
- Extremely comminuted fractures where insufficient holding power of the balloon on the intramedullary canal is probable.
- Delivery sheath is unable to cross fracture site after proper fracture reduction and realignment.

**WARNINGS**

Follow the instructions and warnings provided by suppliers of any cleaning and disinfection agents or equipment used.

High Alkaline solutions can damage instruments with Aluminum components.

Avoid exposure to hypochlorite solutions, which promote corrosion.

Do not reuse or attempt to re-sterilize disposable components. Reuse of single use devices creates a potential risk for patient or user which could lead to injury or death.

**INSTRUMENT STERILITY, CLEANING AND DECONTAMINATION**

Clean instruments prior to initial sterilization and as soon as possible after use. Do not allow blood to dry on the instruments. If cleaning must be delayed, place groups of instruments in a covered container or enzymatic solution to delay drying. Wash all instruments, even if they were not used and did not appear to come in contact with blood or saline solution.

Ultrasonic cleanser can be used with hot water per manufacturer's recommended temperature (usually 90-140 degree F or 32-60 degree C)

900072\_F

and specially formulated detergents. Follow manufacturer's recommendations for proper cleaning solution formulated specifically for ultrasonic cleaners. Be aware that loading pattern, instrument cassettes, water temperature and other external factors may change the effectiveness of the equipment.

Washer-decontamination equipment will wash and decontaminate instruments. Complete removal of soil from crevices and serrations depends on instrument construction, exposure time pressure of delivered solution, and pH of the detergent solution, and thus may require brushing.

Be familiar with equipment manufacturer's use and operation instructions. Be aware that loading, detergent, water temperature and other external factors may change the effectiveness of the equipment.

The effectiveness of subsequent decontamination processes depends on prior removal of gross soil and it may be impaired by dried or coagulated protein. Gross soil should be removed under running water using a mechanical aid, such as a brush with rigid nylon bristles.

Care should be taken to avoid splashing and generating aerosols by holding instruments below the surface of the water in a sink into which water is running and continuously draining. Instruments should not be held under a running tap, as it is likely to result in splashing.

Operatives should wear protective equipment including gloves and goggles. Care should be taken to avoid penetrating or cutting injuries. Particular attention should be taken to remove all debris from all cannulations and obscure holes in the instruments.

It is recommended that the instruments, disassembled as required, be decontaminated using an automatic washer disinfection unit utilizing thermal disinfection. This should preferably be of the ultrasonic or continuous tunnel process type. The cabinet type is an acceptable alternative if a continuous process machine is unavailable.

Typical initial cleaning temperature of the instrument should reach a minimum temperature of 160°F (71°C) for a minimum of 3 minutes, 176°F (80°C) for a minimum of 1 minute or 194°F (90°C) for 1 second. Compatible detergents and rinse aids may be used as recommended by the manufacturer of the washer-disinfectant unit. These detergents and /or rinse aids, however, should be of neutral or near neutral pH.

The recommended sterilization cycle for instruments using PreVac is a minimum of 4 minutes at 270°F (132°C) and 30-32 PSIG or Gravity is a minimum of 15 minutes at 270°F (132°C) and 30-32 PSIG.

After cleaning/disinfecting, the disassembled instruments should be reassembled and placed in their proper locations in the instrument cases.

Surgical Instruments and Instrument cases are susceptible to damage for a variety of reasons including prolonged use, misuse, or rough or improper handling. Care must be taken to avoid compromising their exacting performance. To minimize damage and risk of injury, inspect the instrument case and instruments for damage upon receipt and after each use and cleaning. Incompletely cleaned instruments should be re-cleaned. Instruments in need of repair should be set aside for repair service or returned to IlluminOss.

Instruments returned to IlluminOss or its distributors should be cleaned and sterilized prior to shipment. ANSI/AAMI ST35 Safe Handling and Biological Decontamination of Reusable Medical Devices in Health Care Facilities and in Nonclinical Settings provides guidelines for return, or contact IlluminOss or your distributor for instruction.

Only use an instrument for its intended purposes.

WHEN HANDLING SHARP INSTRUMENTS, USE EXTREME CAUTION TO AVOID INJURY. CONSULT WITH AN INFECTION CONTROL PRACTITIONER TO DEVELOP AND VERIFY SAFETY PROCEDURES APPROPRIATE FOR ALL LEVELS OF DIRECT INSTRUMENT CONTACT.

**RECOMMENDED STORAGE**

Store sterile packaged instruments in a manner that protects them from dust, moisture, insects, vermin and extreme heat and humidity. If desired, use instrument trays to contain instruments that are provided in sets. Double wrap instruments according to local procedures using techniques such as those outlined in ANSI/AAMI ST46-1993.

**GENERAL MAINTENANCE**

Regularly check instrument for wear and damage, such as bending or distortion. Cutting edges should have continuous edge and be free of nicks, jaws and teeth should align and engage properly. Movable parts should move smoothly without excessive play and locking mechanisms should fasten securely. Lubricate moving components with water soluble lubricant between uses according to manufacturers instructions.

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U.S.A.