Description

IlluminOss offers sets with a variety of standard, reusable surgical instruments to assist orthopedic surgeons with the implantation procedure for the IlluminOss Photodynamic Bone Stabilization System. The intended purpose of the IlluminOss Reusable Instruments is for use in conjunction with the IlluminOss implants, for obtaining access to the bone, the preparation of the bone canal, the delivery of the implant, and if needed the removal of the implant. The IlluminOss Reusable Surgical Instruments are used only in conjunction with the IlluminOss Photodynamic Bone Stabilization System, and therefore the indications, contraindications, intended patient population, and intended users that govern their use are those of the implant, the IlluminOss Photodynamic Bone Stabilization System.

These reusable surgical instruments must be cleaned and sterilized before use. The information provided is for the reusable surgical instruments only and does not apply to the Photodynamic Bone Stabilization System catheters (implant), or sterile disposable ancillary devices. Refer to the implant Instructions for Use for the implant indications, contraindications, warnings, precautions and risks (ref: 900535 found on www.illuminoss.com).

Contraindications

- Do not use stainless steel instruments in patients who are allergic to nickel

Warnings

- Before completing an Illuminoss implantation refer to implant instructions for use the implant indications, contraindications, warnings, precautions and risks (ref: 900535 found on www.illuminoss.com).
- Before completing an Illuminoss implantation, refer to the appropriate surgical technique guide for proper surgical technique.
- Stainless steel instruments may contain hazardous substance Cobalt in concentrations up to 0.5% w/w within the stainless steel alloy.
- Incompletely cleaned instruments should be re-cleaned.
- Thoroughly clean instruments prior to initial sterilization.
- The Surgical Instrument Sets and their instruments are supplied non-sterile and must be sterilized before use unless otherwise indicated.
- Clean instruments as soon as possible, after use. Do not allow blood to dry on the instruments.
- Do not stack trays during sterilization.
- Automated washer/disinfector systems are not recommended as the sole cleaning method for surgical instruments. Instruments should be cleaned following the manual or combination manual/automated cleaning procedure.
- All instruments should be placed in their designated instrument trays and the tray wrapped with an approved CSR sterilization prior to sterilization. Effective sterilization is predicated on thorough cleaning and drying processes. Failure to do so may compromise the sterilization process and render the processed instruments unsuitable for clinical use.
- Instruments with dents, burs, signs of corrosion or spotting, or raised surfaces which could cause damage to surgical gloves should be set aside for repair service or returned to IlluminOss.
- Do not reuse or attempt to clean/sterilize single use instruments. These instructions are applicable to the IlluminOss Reusable Surgical Instruments, only.

Precaution

- Stainless steel instruments contain nickel within the alloy.
- Use caution when handling sharp instruments, use extreme caution to avoid injury.
- Do not apply excessive torque when using instruments, excessive torque can result in instrument breakage- excessive torque is defined as > 10 degree bending of instrumentation during usage and using sizes greater than the ID of the canal.
- Reamers are intended to be used to clear only soft tissue from the intramedullary canal and not for reaming cortical bone, reaming of cortical bone may result in instrument breakage.
- Complete proper sizing the canal prior to choosing instrumentation. Utilization of instrumentation larger than the canal may result in applying excessive torque to instrumentation which could result in breakage.
- Do not use flexible instruments if they have a permanent bend or kink in them or are not straight.
- Do not use flexible reamers if there are gaps in consistency along the laser cuts.
- Do not use instrumentation in conjunction with high torque adapters used for power drills.

- Refer to Appendix 1 below for dimensional compatibility of instrumentation.
- Flexible reamers and flexible burrs should not be used if the canal contains bone cement, metallic implants or other surgical instruments. Doing so may result in instrument breakage.

Residual Risks include:

- Infection, including wound complications
- Adverse tissue reaction
- Damage to cutaneous nerves, posterior structures, or injury to arteries or other small blood vessels. Loss of sensation, pain, or range of motion in treated anatomy
- Inability to properly use or remove device. Portion of device may remain in patient.
- Extended surgical and anesthesia time
- Failed fixation, discomfort, possible implant removal / revision using alternative fixation methods
- Soft tissue damage and/or bone loss
- For other risks related to the Photodynamic Bone Stabilization System see IFU 900353.

Instrument Sterility, Cleaning, and Decontamination

The Surgical Instrument Sets and their instruments are supplied non-sterile and must be sterilized before use unless otherwise indicated. Thoroughly clean instruments prior to initial sterilization. When handling sharp instruments, use extreme caution to avoid injury. Consult with an infection control practitioner to develop and verify safety procedures appropriate for all levels of direct instrument contact. Instruments returned to IlluminOss or its distributors should be cleaned and sterilized prior to shipment. ANSI/ AAMI ST35 “Safe Handling and Biological Decontamination of Reusable Medical Devices in Health Care Facilities and in Nonclinical Settings” provides guidelines for return or contact IlluminOss or your distributor for further instruction.

Point of Use:	Clean instruments as soon as possible, after use. Do not allow blood to dry on the instruments. Wash all instruments whether or not they were used or inadvertently came into contact with blood or saline solution. • If cleaning must be delayed, place groups of instruments in a covered container with appropriate detergent or enzymatic solution to delay drying. It is recommended to transport the devices from point of use to processing in either the designated instrument tray or the covered container with solution if used, in a specially designated trolley, or similar, by trained staff.
Cleaning:	Manual Cleaning Instructions: A neutral pH / enzymatic detergent is the recommended detergent for this cleaning process. • Immediately after use in the surgical procedure and instrument disassembly, remove as much visible debris as possible from each instrument and its components using moist clean gauze pads or wipes. Instruments should be soaked immediately after use. If the instruments cannot be soaked immediately, they should be wrapped in a clean moist towel. Soiled instruments must be kept moist to prevent the soil from drying. • Prepare an enzymatic cleaning solution with warm clean utility water (35-43°C) per the instructions of the solution manufacturer. Perform an enzymatic soak for a minimum of 20 minutes making sure the instruments are completely immersed in the solution. Actuate all moving mechanisms a minimum of 5 times (5x) to ensure the enzymatic solution contacts all parts. • While in the soak solution, use a soft brush or clean cloth to manually remove all the exterior soil. Use appropriately sized brushes to thoroughly clean the entire length of any lumens including all other challenging design features such as holes, hinge/mated surfaces, crevices, serrations, etc. Scrub each challenging design feature a minimum of 1 minute, or approximately 5 times (5x). • After the enzymatic soak, rinse instruments thoroughly with clean warm utility water (35-43°C), taking care to flush all lumens and crevices. Rinse thoroughly for at least one minute under a continuous stream of clean water ensuring all surfaces are rinsed and the instruments visibly clean. If rinse water does not run clear or there are any visual signs of soil, repeat cleaning process. • Prepare a fresh enzymatic cleaning solution with warm clean utility water (35-43°C) per the instructions of the solution manufacturer. This solution will be used in the ultrasonic bath. The instruments should be fully immersed in the cleaning solution and all mechanisms actuated a minimum of 5 times (5x). Run the ultrasonic cleaning unit for a minimum of 15 minutes. • Place the instruments in clean warm (35-43°C) critical water and again scrub the instruments, ensuring all the surfaces are scrubbed, including the lumens, interfaces, and crevices. Scrub each challenging design feature a minimum of 1 minute, or approximately 5 times (5x). Rinse each instrument again for at least one minute under clean warm (35-43°C) running critical water (at a minimum flow rate of 1 Liter/minute), including all design

	features. Ensure that all surfaces are free of visible contamination or any traces of the cleaning solution (e.g. foam) and the fluid runs clear. • Dry the instruments with a clean soft cloth and use clean instrument air to dry the lumens, holes, and interfaces. • Perform a visual inspection on the instruments and verify that they are clean, dry, and in proper working order. If the instruments are not visually clean, repeat the manual cleaning process. Automated Cleaning Instructions (Optional): • Automated washer/disinfector systems are not recommended as the sole cleaning method for surgical instruments. Instruments should be cleaned following the manual or combination manual/automated cleaning procedure. • Devices should be thoroughly inspected prior to sterilization to ensure effective cleaning. <table border="1"> <thead> <tr> <th>Step</th><th>Description</th></tr> </thead> <tbody> <tr> <td>1</td><td>5 min pre-wash with cold tap water</td></tr> <tr> <td>2</td><td>10 min detergent wash with a neutral pH / enzymatic detergent at 60°C (an alternate that may be used is an alkaline pH / non-enzymatic detergent)</td></tr> <tr> <td>3</td><td>1 min rinse with hot utility water</td></tr> <tr> <td>4</td><td>Thermal disinfection $A_0 = 3000$ at 93°C (approx. 10 min)</td></tr> <tr> <td>5</td><td>20 min hot air drying at 110°C</td></tr> </tbody> </table> Upon completion, unload the washer/disinfector. Perform a visual inspection of the instruments and verify that they are clean, dry, and in proper working order. If the instruments are not visually clean, repeat both the manual and automated cleaning process. Note: The washer/disinfector manufacturer's instructions should be strictly adhered to. Use only cleaning agents recommended for the specific type of automated washer/disinfector. A washer/disinfector with approved efficacy (e.g. CE mark, FDA approval, and validation according to ISO 15883) should be used.	Step	Description	1	5 min pre-wash with cold tap water	2	10 min detergent wash with a neutral pH / enzymatic detergent at 60°C (an alternate that may be used is an alkaline pH / non-enzymatic detergent)	3	1 min rinse with hot utility water	4	Thermal disinfection $A_0 = 3000$ at 93°C (approx. 10 min)	5	20 min hot air drying at 110°C
Step	Description												
1	5 min pre-wash with cold tap water												
2	10 min detergent wash with a neutral pH / enzymatic detergent at 60°C (an alternate that may be used is an alkaline pH / non-enzymatic detergent)												
3	1 min rinse with hot utility water												
4	Thermal disinfection $A_0 = 3000$ at 93°C (approx. 10 min)												
5	20 min hot air drying at 110°C												
Sterilization:	• All instruments should be placed in their designated instrument trays and the tray wrapped with an approved CSR sterilization prior to sterilization. Effective sterilization is predicated on thorough cleaning and drying processes. Failure to do so may compromise the sterilization process and render the processed instruments unsuitable for clinical use. • Wrap Instrument Case in an approved CSR sterilization wrap according to local procedures using techniques such as those outlined in AAMI ST79. Sterilize in an approved sterilizer. • Do not stack trays during sterilization. • Sterilize Instruments using PreVac steam autoclave or Gravity steam autoclave: 1. For PreVac steam autoclave, sterilize instruments for a minimum of 4 minutes at 270°F (132°C) and at cycle completion, dry for a minimum of 20 minutes at 270°F (132°C). OR 2. For PreVac steam autoclave, sterilize instruments for a minimum of 3 minutes at 134°C (273°F) and at cycle completion, dry for a minimum of 16 minutes at 134°C (273°F). OR 3. For Gravity steam autoclave, sterilize instruments for a minimum of 15 minutes at 270°F (132°C) and at cycle completion, dry for a minimum of 15 minutes at 270°F (132°C).												
After Sterilization	• After removal from the sterilizer, move to an area where the trays are allowed to cool. When cool, place in a sterile store until delivery to the point of use. The storage area needs to be a restricted area and away from windows or traffic. The sterilized trays should be handled as little as possible. • Sterile trays should be stored in a clean, dry and well-ventilated area with open racks for air circulation and moderate temperature and humidity ($\leq 60\%$ RH). Avoid storing in closed cupboards or on the floor. Sterile trays should be stored away from direct sunlight and contamination with dust, moisture or gases, e.g., exhaust fumes. They should be placed securely away from traffic to prevent accidental physical damage and stored according to local procedures. If the tray is dropped after sterilization, repeat inspection and sterilization. Do not store sterile trays longer than recommended by the CSR sterilization wrap used or local procedures. • It is recommended to transport the sterile devices between the decontamination unit and point of use/storage in a specially designated clean trolley, or similar, by trained staff. • Immediately before use, inspect the sterile wrapping used on the tray. If any holes, tears, stains, or damage is observed, repeat device sterilization before use.												
Inspection:	• Inspect the instrument case and instruments for wear or damage upon receipt and after each use and cleaning. Cutting edges should have continuous edge and be free of nicks. Jaws and teeth should align and engage properly. Movable parts should move smoothly without excessive play and locking mechanisms should fasten securely. Do not use worn or damaged instruments. Worn or damaged instruments, including flexible instruments with bends, kinks, or gaps in consistency along the laser cut, could indicate the instrument is at the												

	end of its life and should not be reprocessed and used. Contact IlluminOss immediately for return, repair, or replacement of damaged or worn instruments. • Instruments with dents, burs, signs of corrosion or spotting, or raised surfaces which could cause damage to surgical gloves should be set aside for repair service or returned to IlluminOss. • Incompletely cleaned instruments should be re-cleaned. • Do not reuse or attempt to clean/sterilize single use instruments. These instructions are applicable to the IlluminOss Reusable Surgical Instruments, only. • Any serious incident that has occurred in relation to the device should be reported to the manufacturer via www.illuminoss.com and to the competent authority of the Member State in which the user is established. Reporting to the Authorized Representative is one means of doing so.
--	--

MD**REF**

SY-3000 SY-4000 SY-5000 SY-9100

CT-0239 CT-0241 CT-0242 CT-0186 CT-0188 CT-0196 05-0901 CT-0129 CT-0131 CT-0238 CT-0278
 CT-0235 CT-0268 CT-0037 CT-0038 CT-0039 CT-0040 CT-0041 CT-0042 CT-0151 CT-0146 CT-0143 CT-0161
 CT-0156 CT-0153 CT-0162 CT-0163 CT-0199 CT-0127 CT-0195 CT-0197 CT-0301 CT-0267 CT-0266 CT-0206
 CT-0207 CT-0231 CT-0232 CT-0233 CT-0225 CT-0226 CT-0227 CT-0228 CT-0229 CT-0237 CT-0270 CT-0216
 CT-0305 CT-0306

ILLUMINOSS MEDICAL INC.
 993 Waterman Avenue
 East Providence, RI 02914
 U.S.A.
 +1 (401) 714-0008

EMERGO EUROPE
 Westervoortsedijk 60,
 6827 AT Arnhem
 The Netherlands

**2797**

1

Stainless steel instruments may contain hazardous substance Cobalt in concentrations up to 0.5% w/w within the stainless steel alloy. Cobalt is classified as a category 1B carcinogen and category 1B reproductive toxin, meaning it has presumed carcinogenicity and reproductive toxicity. Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or reproductive effects.

The typical maximum speed of common cordless surgical rotary handpieces (power drill) ranges up to about 1500 RPM. This occurs when in “drill mode” or using drill attachments, but this speed can be reduced to about 270 RPM if a “ream mode” switch option is provided with the handpiece or when using reamer attachments. The RPM of the reamer is not critical to the clearing of soft tissue from the intramedullary canal in preparation for placement of the implant. IlluminOss recommends the use of the flexible reamers in drill mode as the instrument is intended to clear soft tissue from the intramedullary canal to open the space for the implant. The flexible reamers are not intended to cut cortical bone in the intramedullary canal.

Instrument Product Name	Product Code / Catalog #	Instrument End Connector Type	Power Handpiece Attachment Type(s) Required
4.5mm Front Cutting Flexible Burr	CT-0239	Hex (6 flat) drive, 5.5mm	1/4" (6.4mm) Chuck or larger. Keyed or keyless chucks available. Chuck aka "Jacobs chuck".
5.5mm Side Cutting Flexible Burr	CT-0241		
6.5mm Side Cutting Flexible Burr	CT-0242		
AO Adapter	CT-0238	Snap-lock shank for AO small drill attachment	AO, small drill <i>Accepts Synthes tools with the appropriate snap-lock shank</i>

Drill Bit, 8mm	CT-0216	3-flat drive, 7.0mm	3/8" (9.5mm) Chuck or larger. Keyed or keyless chucks available. Chuck aka "Jacobs chuck"
6.0mm Flex Reamer	CT-0225	Hudson / Modified Trinkle (quick disconnect) Connector design includes distal 7.3mm 3-flat portion	Hudson / Modified Trinkle (accepts accessories such as reamers, drill bits, and automatic screwdrivers with Hudson or modified Trinkle fittings; industry standard design) 3/8" chuck also designed to secure connector which has a 7.3mm, 3-flat design.
6.5mm Flex Reamer	CT-0226		
7.0mm Flex Reamer	CT-0227		
7.5mm Flex Reamer	CT-0228		
8.0mm Flex Reamer	CT-0229		
Head Reamer (18mm)	CT-0231	Hex (6 flat) drive, 5.5mm	1/4" (6.4mm) Chuck or larger. Keyed or keyless chucks available. Chuck aka "Jacobs chuck".
Head Reamer (20mm)	CT-0232		
Head Reamer (22mm)	CT-0233		
Cannulated Drill Bit, 5mm	CT-0037	3-flat drive, 5.0mm	1/4" (6.4mm) Chuck or larger. Keyed or keyless chucks available. Chuck aka "Jacobs chuck".
Cannulated Drill Bit, 6mm	CT-0038	3-flat drive, 6.0mm	
Cannulated Drill Bit, 7mm	CT-0039	3-flat drive, 7.0mm	3/8" (9.5mm) Chuck or larger. Keyed or keyless chucks available. Chuck aka "Jacobs chuck"
Cannulated Drill Bit, 8mm	CT-0040	3-flat drive, 6.3mm	1/4" (6.4mm) Chuck or larger. Keyed or keyless chucks available. Chuck aka "Jacobs chuck".
Cannulated Drill Bit, 9mm	CT-0041		
Cannulated Drill Bit, 10mm	CT-0042		
2.0mm Threaded Tip Guidewire, 18"	CT-0143	2.0mm diameter wire.	1/4" (6.4mm) Chuck. Keyed or keyless chucks available. Chuck aka "Jacobs chuck". Wire/Pin Collets.
2.0mm Threaded Tip Guidewire, 12"	CT-0146		
2.0mm Threaded Tip Guidewire, 6"	CT-0151		
2.5mm Threaded Tip Guidewire, 18"	CT-0153	2.5mm diameter wire.	
2.5mm Threaded Tip Guidewire, 12"	CT-0156		
2.5mm Threaded Tip Guidewire, 6"	CT-0161		

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